

# BALUM

# 61 Times

# REUSABLE

## BASKA MASK


Manufactured and Sponsored by  
**LOGIKAL HEALTH PRODUCTS Pty Ltd**  
UNIT 3, 18-20 ACCOLADE AVENUE  
MORISSET NSW 2264 AUSTRALIA

INTERNATIONALLY PATENTED  
AND PATENTS PENDING


[www.baskamask.com](http://www.baskamask.com)

## INSTRUCTIONS FOR USE

To ensure integrity of mask, with gloved hands, occlude the airway opening of the connector end of the device with one thumb, hold the mask head with the other hand and place the thumb over the airway opening of the mask to seal. Apply pressure for 5 seconds to confirm there is no leak in the device.




**⚠ Do not use if leak is detected ⚠**



Lubricate the entire mask part of the device liberally with a water soluble lubricant. This is the key to a successful insertion and placement.

Mildly lubricate both the upper and lower lips of the patient. Place the head and neck in slightly extended position. With the mouth open insert the device all the way down holding it at it's stem. Do not pull on the lower jaw. Allow the lower lip to cover over the teeth during insertion. Do not twist or turn the device while inserting. If required, pull the tab gently to help negotiate the palato-pharyngeal curve. Release the tab as soon as the tip of the mask has passed around the curve.

Ensure the mask has been advanced sufficiently until resistance is encountered. In this position the tip of the mask is engaged into the upper end of the oesophagus. Check that the airway is clear, by observation if the patient is breathing spontaneously or by gently ventilating if the patient is apnoeic. It may occasionally be necessary to withdraw the mask very slowly, 3-5mm at a time until a clear airway is obtained. At this point secure the mask in position.



Complete and check the seal by closing the breathing circuit spill valve and applying pressure to the reservoir bag.

If regurgitation is anticipated or is a high risk, suction should be connected to the suction port of the device before insertion and be left running continuously until the mask has been positioned correctly. Suction should also run continuously during removal of the device at the end of the procedure. During the procedure suction can be applied intermittently as required to clear any fluids or secretions.

## CLEANING and STERILISATION:

Refer to BALUM Reusable BASKA Mask Record Card supplied along with the product in the pouch.

CE  
0805



## CAUTION:

**Only to be used by an adequately trained person.**  
**Do not use if the package or product has been damaged.**  
**Do not use if the product is distorted in shape.**  
**Do not store at extreme temperature, humidity or UV conditions.**

EUROPEAN AUTHORISED REPRESENTATIVE

**PRIORY ANALYSTS Ltd**

THE PINNACLE  
160 MIDSUMMER BOULEVARD  
MILTON KEYNES, MK9 1FF  
UNITED KINGDOM

**LATEX**

## REUSABLE – STERILISE BEFORE USE

# 61